



## Cavity Ball / Bola para cavidades




D 001-008 D 001-010 D 001-012 D 001-014 D 001-016 D 001-018 D 001-021 D 001-024 D 001-027 D 001-033 D 001-035

## Cavity Ball / Bola para cavidades




D 002-008 D 002-010 D 002-012 D 002-014 D 002-016

## Long Cavity Ball / Bola larga para cavidades



D 001L-012 D 001L-014 D 001L-016 D 001L-018 D 001L-023

## Inverted Cone / Cono invertido




D 010-010 D 010-016 D 019-010 D 019-012 D 019-014 D 225-016 D 225-018

## Barrel / Barril




D 038-037

## Occlusion / Oclusión




D 507-023 D 466-023

## Cavity Pear / Pera para cavidades



D 233-009 D 233-010 D 234-012 D 234-014 D 234-018

## Micro Prep 80 micron / Micropreparación 80 µm



M 276-009 M 276-012 M 276-014 M 536-007 M 137-007

## Fissure / Fisuras




D 465-006

## Interdental / Interdental




D 465-016

## Short Shank / Vástago corto (Bambino)




M 001-010 M 001-014 M 156-012 M 196-012 M 233-012 M 288-012

## Safety End / Punta de seguridad




F 141-014 F 142-016 F 298-018 F 299-021

## Wheel / Rueda




D 040-018 D 068-040 D 068-052

## Separator / Separador



D 199-012s D 167-010s

## Rounded Edge Cylinder / Cilindro canto redondeado




D 156-008 D 156-010 D 156-012 D 156-014 D 157-010 D 157-012 D 157-014 D 158-010 D 158-012 D 158-014 D 158-018

## Long Cylinder / Cilindro largo




D 158L-014

## End Cutting Only / Corte sólo punta




D 150-010 D 150-012

## Round End Cylinder / Cilindro punta redonda



D 137-008 D 137-010 D 137-014 D 140-012 D 140-014 D 141-010 D 141-012 D 141-014 D 142-012 D 142-014 D 143-014

## Amalgam Reducer / Reductores de amalgama



Zyl 5 Zyl 7 Zyl 9 Kon 8 Kon 10


■ normal grit / grano normal  
■ coarse grit / grano grueso

■ fine grit / grano fino  
■ super fine grit / grano superfino


each pack contains 5 instruments / estuches de 5 instrumentos

scale 1:1 / escala 1:1


## Flame / Llama




## Pointed Torpedo / Torpedo en punta



## Torpedo / Torpedo




## Tapered Torpedo / Torpedo cónico



## Flat Taper / Cono punta plana




## Round End Taper / Cono punta redonda




## Safe End Taper (endo) / Cono punta inactiva (endo)




## Zirconium Bur / Fresa para circonio




## Flat Taper, Rounded End / Cono canto redondeado




## Egg / Ovoide



## Pointed Egg / Ovoide en punta



## Bud / Granada



■ normal grit / grano normal  
■ coarse grit / grano grueso

■ fine grit / grano fino  
■ super fine grit / grano superfino

each pack contains 5 instruments / estuches de 5 instrumentos


scale 1:1 / escala 1:1

# upper-mant *Supra*


|                                                                                                                                                                                      |                                                             |                                                                                                      |                                                    |                                                                                                                                 |  |                                                               |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------|------------------------------------------------------------------------------------------------------|----------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------|--|---------------------------------------------------------------|
| <b>Ball / Bola</b><br><br>S 001-014 S 001-021<br>                                                                                                                                    | <b>Pear / Pera</b><br><br>S 234-012 S 234-014 S 234-018<br> | <b>Bud / Granada</b><br><br>S 274-020<br>                                                            | <b>Egg / Ovoide</b><br><br>S 277-020 S 277-023<br> | <b>Cylinder Rounded Edge / Cilindro canto redondeado</b><br><br>S 156-012 S 156-014 S 157-012 S 157-014 S 158-012 S 158-014<br> |  |                                                               |
| <b>Torpedo / Torpedo</b><br><br>S 288-012 S 288-014 S 289-010 S 289-012 S 289-014 S 289-016 S 290-012 S 290-014 S 290-016<br>                                                        |                                                             |                                                                                                      |                                                    | <b>Round End Cylinder / Cilindro punta redonda</b><br><br>S 140-012 S 140-014 S 141-012 S 141-014<br>                           |  |                                                               |
| <b>Tapered Torpedo / Torpedo cónico</b><br><br>S 297-012 S 297-014 S 297-016 S 298-014 S 298-016 S 298-018 S 298-020 S 298-022 S 299-012 S 299-014 S 299-016 S 299-018 S 299-022<br> |                                                             |                                                                                                      |                                                    |                                                                                                                                 |  | <b>Flame / Llama</b><br><br>S 250-012 S 250-014 S 250-016<br> |
| <b>Flat End Taper / Cono punta recta</b><br><br>S 171-016 S 171-023 S 172-016<br>                                                                                                    |                                                             | <b>Rounded Edge Taper / Cono canto redondeado</b><br><br>S 544-016 S 545-016 S 545-023 S 546-016<br> |                                                    | <b>Round End Taper / Cono punta redonda</b><br><br>S 198-012 S 198-018 S 199-012 S 199-018<br>                                  |  |                                                               |

## Double Layer Diamond Burs Fresas de Doble Capa de Diamante

### conventional technology



### SUPRA



### SUPRA Diamonds have a Longer Cutting Power

Cutting Capacity of a Conventional Diamond Bur

SUPRA Diamond cuts for a much longer time



■ normal grit / grano normal  
■ coarse grit / grano grueso

■ fine grit / grano fino  
■ super fine grit / grano superfino


each pack contains 5 instruments / estuches de 5 instrumentos

scale 1:1 / escala 1:1



UPPER-MAXI S.L.

## Finishers / Acabadores




## Endo Bur / Fresa Endo




# TUNGSTEN-CARBIDE BURS FRESAS DE CARBURO-TUNGSTENO


## Round Carbide Bur / Bola de tungsteno




## Crown Cutter / Cortador de coronas




## Orthodontic Bur / Ortodoncia




## Arkansas / Arkansas




## Composite & Amalgam Polishers / Pulidores de composite y amalgama




## Prophylaxis / Profilaxis



## Ceramic Polishers / Pulidores para cerámica



## Laboratory Tungsten-Carbide Burs / Fresas de carburo-tungsteno para laboratorio





Group D: For plaster and acrylic.  
Grupo D: Para yeso y resina.


Group F: For acrylic, metal and noble metal.  
Grupo F: Para resina, metal y metal noble.

Group G: For noble metal, acrylic and plaster.  
Grupo G: Para metal noble, resina y yeso.

## Laboratory Polishers / Pulidores de laboratorio



## Toledo Elastic Diamond Disks / Discos Toledo flexibles de diamante







FG : friction grip / turbina

RA : contra-angle / contra-ángulo

HP: hand piece / pieza de mano

scale 1:1 / escala 1:1

## SURGICAL BURS FRESAS DE CIRUGIA



Each pack contains 3 instruments, except implantology drills.  
Cada estuche contiene 3 instrumentos, excepto fresas de implantología.

FG : friction grip / turbina

RA : contra-angle / contra-ángulo

HP : hand piece / pieza de mano